

Jerry Zion presents:

Medical Device Quality Assurance Testing: Best Practices For Patient Risk Reduction

Sponsored by: Fluke Biomedical





Wednesday, June 28, 2017, 2:00pm ET

PRESENTER: JERRY ZION



Jerry Zion, Global Training Manager at Fluke Biomedical, is a renowned veteran in the medical device industry, and has more than 35 years of experience working in many capacities for hospitals and various medical device manufacturers. Jerry has spent the last decade at Fluke Biomedical as marketing manager, product manager, and pioneering the global training program. Jerry is an AAMI certified biomedical equipment technician, and holds a Bachelor of Science degree in Electrical Engineering Technology from Purdue University. Jerry also obtained a Master's in Management in Science and Technology from the Oregon Graduate Institute.

WEBINAR AGENDA: During this 60-minute presentation join Jerry Zion, Global Training Manager at Fluke Biomedical, as he discusses applying the Hypocratic Oath: "first, do no harm", Reducing patient risk using an effective quality management program; Sources of patient risk; Ensuring the clinical staff know how to get the most from their medical devices; Maintaining the medical device history record and compliance to GMP (repairs, inspections, etc.) as well as Keeping instrumentation current with Innovations in medical devices.

At the end of the session, participants will:

- 1) Understand the sources of patient risk.
- 2) Understand that human error is not managed using instrumentation or technology alone.
- 3) Understand how learning about the ways others have managed patient risk sources instructs/improves our own plans.
- 4) Understand the importance of maintaining the medical device history record is a responsibility the hospital Biomedical/Clinical Engineering Department cannot avoid under the USA FDA GMP 21 CFR requirements, even when repair, calibration services are outsourced.
- 5) Understand how establishing meaningful key performance indicators (kpis) and keeping daily visual management up to date helps us and visitors understand our role better, and our value in patient risk reduction.

Webinar Wednesday would like to thank our sponsor Fluke Biomedical!



Fluke Biomedical is the premier, global provider of test and measurement equipment and services to the healthcare industry. We serve biomedical engineers, quality-assurance technicians, medical physicists, oncologists, radiation-safety professionals and are continually expanding our range of solutions to a broader range of health and safety professionals.

Learn more at www.flukebiomedical.com



DON'T MISS THE WEBINAR WEDNESDAY'S Q&A

Our presenters will end today's webinar with a Q&A session. If you have a question, please submit it early by emailing webinar@mdpublishing.com.

All questions will be addressed at the end of today's presentation or the presenter will follow up with you offline.



SMILE - SNAP - AND SHARE!

We love seeing photos of the Healthcare Community! Send a photo of you or your department participating in this week's Webinar Wednesday to webinar@mdpublishing.com. Submissions will be used for the magazine, website and social media so smile, snap a photo and share with us!

Use our hashtag

#WEBINARWEDNESDAY

on Facebook and Twitter

"I find Webinar Wednesday's highly educating and informative; I simply love it. It makes me a smarter BioMed!"
-Alero O.

"Webinar Wednesday's are an excellent resource for continuing education credits. We have found them to answer some of the questions we have had for specific equipment management, and the discussions on regulations and latest rules have been useful as well. Keep up the good work!"

- Leroy S.



Medical device quality assurance

Best practices for patient risk reduction

The Hypocratic Oath



----first, do no harm----

The simple act of connecting a patient to a medical device exposes that patient to potential harm if the medical device fails. This is a risk that must be managed

How do we do that?



Quality management program



- Ensure the ongoing quality, safety and effectiveness of medical devices
- Avoid medical device failures during critical medical procedures
- Improve clinical effectiveness
- Reduce total cost of device ownership
- Uphold staff morale and professionalism through positive patient experience
- Reduce patient risk and improve patient confidence
- Comply with regulations and OEM recommendations







Medical device management plan: clinical factors



Device function

- What function does the equipment perform in a clinical environment?
 - Highest risk → life-support devices
 - Lower risk → with non-invasive, diagnostic devices

Risk of misuse or failure

- What are the possible consequences to the patient of a device malfunction or misapplication?
 - Range from "no significant risk" to death

Mission criticality

- What is the impact on overall hospital patient care or patient flow when the medical device in question is not available?
 - Critical, important, non-critical





Patient risk reduction





How do we ensure the continuous safety and effectiveness of medical devices in an systematic, repeatable way?


Medical device technology should:

- Help sustain high quality of care provided to each patient
- Facilitate faster/easier diagnoses to improve patient care/outcomes
- Be safe to use at all times for both patients and staff
- Maintain functionality and effectiveness throughout the life of the device



Quality assurance: maintenance factors





Manufacturer maintenance requirements

- Recommendations based on device type, design, and the components inside
- Compliance with standards

Equipment maintenance history

- How prone to failure is this device or group of devices?
- Maintenance sensitivity: device failures more likely in the absence of scheduled maintenance and testing




Organization and keeping track



- Computerized maintenance management systems (CMMS)
 - Many providers worldwide
 - Some in-house-designed (Microsoft Access, etc.)
 - Keep track of medical device inventory
 - Help schedule and assign inspection work
 - Keep medical device history (repairs, inspections, etc.)
- Risk management databases
 - ECRI Institute Health Devices Alerts
 - ASHE Inspection system


All USA hospitals have Clinical Engineering services and most use a CMMS to track/organize/analyze and plan with quality assurance data






Accuracy vs Precision






What is Measurement Error?





What is Measurement Error?






What is Uncertainty?




- A number that characterizes the dispersion of values attributed to a measurand, based on the information at hand
- An interval about the measured value that we are reasonably confident that the true value lies within
- How confident? Depends on what level you need
 - For most Calibration Labs, it is expressed at about 95%



Measurement Quality





- Measurement Decision Risk
- If your uncertainty is large compared to the specified requirement, how confident are you in your declaration of In Tolerance or Out of Tolerance?



What is Traceability?



Traceability noun | tra-cea-bil-ity | \trā-sə-'bi-lə-tē\

Property of a measurement result whereby the result can be related to a reference through a documented unbroken chain of calibrations, each contributing to the measurement uncertainty.

JCGM 200:201 International vocabulary of metrology – Basic and general concepts associated terms (VIM) 3RD edition

©2014 Fluke Biomedical www.flukebiomedical.com



What is Traceability?



- Step-by-step transfer process
- Multiple steps extending to the SI
- Each transfer increases the uncertainty

©2014 Fluke Biomedical www.flukebiomedical.com



Traceability – Unbroken Chain




Josephson Junction Array at Fluke, EVT USA



0.02 μV/V **DVMF**









0.1%



Summary



- Metrology underpins our technology and affects patient risk
- Maintain traceability
 - Measurements not traceable? They are not VALID!
- Measurements matter
- Don't know where to start?
 - See bipm.org
 - See FlukeBiomedical.com

©2014 Fluke Biomedical www.flukebiomedical.com



Keeping up with the pace of change



- Innovations in medical devices change the way independent evidence of safety and effectiveness for clinical use is assessed
 - Example: CT scanners with increasing number of "slices" to improve resolution of images and accuracy of diagnosis
- Test instruments and systems change to accommodate innovation and new functionality, streamline data collection and analysis
 - Example: Defibrillator testing
 - Innovation: Impedance sensing functionality to ensure the energy selected is the energy delivered
 - Required change: Use of multiple test loads during testing, including the IEC/AAMI/ANSI 50 ohm test load





Assessing effectiveness



- Quantifiable outcomes:
 - Benchmarked schedule of preventive maintenance results in completion of X% of the plan (e.g., 96% PMs completed)
 - Medical device up-time/availability for use improves by X% (e.g., 25% improvement in up-time)
 - Cost of repair decreased by X% (e.g., 30% reduction in cost of repair)
 - "Spare" inventory of medical devices reduced by X% because main-line medical device inventory is more available (e.g., 30% reduction in medical device "spares" due to improvement in up-time)
- Key performance indicators (KPI's) of performance/benefit of the medical device maintenance team
- CMMS database helps capture/analyze these metrics and helps to report quantified measures of success against plan and/or benchmark

Benchmarking





Growing pressures to reduce cost through efficiency and effectiveness

Medical device quality assurance can quantify, measure and drive efficiency and effectiveness

3

Benchmarking helps determine progress

- •Prove efficiency and effectiveness while still allowing for and evidence-based approach
- •Show the value delivered by investment in these efforts



AAMI Benchmarking resources



AAMI's Benchmarking Solution

Developed by clinical engineering experts, ABS provides detailed benchmarking information on all the important topics:

- Budgeting levels for maintenance, office expense, training, computers/software, and more
- Staffing levels allocated to management, repair, and other responsibilities
- Cost of maintenance contracts
- Hours spent on maintenance and repair
- Personnel qualifications and responsibilities
- Reporting structures
- Number of devices maintained by a CE program
- Percentage of scheduled inspections that identify a need for corrective maintenance

<u>Designing and Developing CE</u>
<u>Department</u> *BI&T* November/ December 2011

How to Use Financial Benchmarks BI&T September/October 2011

Staffing Metrics: A Case Study BI&T July/August 2011

Benchmarking Basics: Cautions and Precautions BI&T May/June 2011

Benchmarking for More, 24x7 Magazine April 2011

AAMI's Benchmarking Solution
Enhanced with New Features AAMI
News, Sept. 2010

Analysis of Cost of Service Ratio and Other Metrics BI&T, July/Aug 2010



It's up to you---



- How well do you understand the goals and objectives of your MDQA program?
- How do you want the medical device to work if you, or someone you love is the next patient to be connected to it?
- Evidence is still a preeminent guide to determine test frequency and scope (*where allowed)





Opportunity for global collaboration



The growing global dedication to medical device quality assurance represents opportunity for global leadership through mentoring

- Medical device safety and performance standards (regulation)
- Clinical Engineering departments and their medical device management plans and experiences
- Interest and partnership in medical device quality assessment and reporting in other countries





Questions?



For more information on radiant warmer or infant incubator testing, contact Fluke Biomedical today:

- Email: sales@flukebiomedical.com

-Toll free: (800) 850-4608

-Direct: (440) 248-9300

-Fax: (440) 349-2307



Want to learn more?





Advantage Training: a courtesy service from Fluke Biomedical

Check out our on-demand training in the Advantage Training Center for more helpful tutorials on this and other subjects.

ATTEND **MDEXPO** FOR FREE!

USE THE VIP PASS BELOW, COMPLIMENTS OF



Biomedical

BOOTH #601

HOTEL BUENA VISTA • ORLANDO, FL • OCTOBER 5-7, 2017

✓ EDUCATION

✓ NETWORKING

√ KEYNOTE

✓ EXHIBIT HALL









REGISTER FOR FREE WITH PROMO CODE:

[**17MDE601**]

View schedule, classes and register online at

WWW.MDEXPOSHOW.COM/ORLANDO

"The MD Expo was excellent. The content was on point and relevant for the attendees. The networking opportunities were excellent." - Christopher N., Corporate Director

ADVANCING THE BIOMEDICAL/HTM PROFESSIONAL

*ONLY valid to those individuals employed within a hospital, health care facility, or who are active military or students. **Proof of employment required.**

Shop multiple vendors in a single online marketplace.

Easily search 1,000s of listings...



New equipment added daily!



THREE-HEAD LED PAX LIGHT GOVERNMENT LIQUIDATION



PATIENT WARMING SYSTEM ZRG



RITTER MIDMARK EXAM TABLE EQUIP STAT



AESPIRE ANESTHESIA MACHINE
MEDIPRO INTERNATION LLC



GE LIGHTSPEED PRO 16 CT INTERMED GROUP



PHILIPS ICT 256 SLICE MTCMED



VENDORS: \
START A FREE
90-DAY TRIAL
TODAY!

MEDICALDEALER.COM

The Simple Way To Buy, Sell And Review



| | | <u> </u> |
|---|--|----------------------------------|
| | | |
| | AULIUS AND THE CASE OF THE CAS | Caraca Adams Constitution at the |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| ; | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| ر و روی چار و رایان با داده با است کالی این فیلستان و در سال و بازنشان و در سالان | | OF SATE VALATO VALATO |
| | | |
| | | |
| | | |
| | | |
| | | |
| × | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |